

# OXX 10

## Adult nasal cannula



Disposable nasal cannula for oxygen therapy.

### Specific Characteristics

Helps minimize irritation and soreness  
 Designed for long-term oxygen patients with curved, non-flared nasal prongs

### Technical data

Length	2 ± 0,1 m
Material	PVC

Class IIa Med.Device according to Dir. 93/42/CEE

OXYGEN NASAL CANNULA  
 OX20110A



CND Classification	R03010203
Registration number	361047

Distributed by Spencer Italia s.r.l.

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)