

OXX 30

10 Anticrushing oxygen tubes



Disposable anti-crushing tube for oxygen therapy.

Specific Characteristics

Anti-crush structure
Easy to connect to a hose holder and to the oxygen mask
The kit includes 10 tubes

Latex Free

Technical data

Lenght	2m
Material	PVC

Class IIa Med.Device according to Dir. 93/42/CEE

10 ANTI-CRUSHING OXYGEN TUBES

OX20131A



CND Classification	R03010204
Registration number	199307

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it