

OXYMASK 150

High concentration oxygen masks



Adult disposable high concentration face mask for oxygen therapy.

Specific Characteristics

Adjustable elastic
Semi transparent shell to improve patient monitoring
Anti-crush tube
Flexible aluminum nose clip to increase stability
Comes with reservoir

Latex Free

Technical data

Material	PVC
Tube length	2 m
N° of masks	10 pcs

Class IIa Med.Device according to Dir. 93/42/CEE

OXYMASK 150 HIGH
CONCENTRATION ADULT, 10 PCS
OX20134A



CND Classification	R03010201
Registration number	104801

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it