

OXYMASK 50

Kit 10 pediatric medium concentration oxygen masks



Disposable face masks for oxygen therapy.

Specific Characteristics

Adjustable elastic
Semi transparent shell to improve patient monitoring
Anti-crush tube
Flexible aluminum nose clip to increase stability

Latex Free

Technical data

Dimensions	100x85x55 ± 10 mm
Material	PVC
Weight ^(1mask + tube)	70 g
Tube length	2 m

Class IIa Med.Device according to Dir. 93/42/CEE

OXYMASK 50 MEDIUM
CONCENTRATION PAED., 10 PCS
OX20161A



CND Classification	R03010201
Registration number	976287

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it