

OXYMASK 100

Kit 10 Variable concentration oxygen masks



Masks for oxygen therapy with venturi system for variation of oxygen concentration.

Specific features

Adjustable elastic
Semi transparent shell to improve patient monitoring
Anti-crush tube
Flexible aluminum nose clip to increase stability

Latex Free

Technical data

Materials	PVC/PP
Tube length	2 m
Corrugated Tube length	15 cm
Obtainable oxygen concentrations	24/28/31/35/40/50/60%

Class IIa Med.Device according to Dir. 93/42/CEE

OXYMASK 100 VARIABLE
CONCENTRATION ADULT 10PCS
OX20170A



Rev.1 (15/01/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it