

## OXYMASK 100

### Kit 10 pediatric Variable concentration oxygen masks



Masks for oxygen therapy with venturi system for variation of oxygen concentration.

#### Specific features

Adjustable elastic  
Semi transparent shell to improve patient monitoring  
Anti-crush tube  
Flexible aluminum nose clip to increase stability

**Latex Free**

#### Technical data

Materials	PVC/PP
Tube length	2 m
Corrugated Tube length	15 cm
Obtainable oxygen concentrations	24/28/31/35/40/50/60%

Class IIa Med.Device according to Dir. 93/42/CEE

OXYMASK 100 VARIABLE  
CONCENTRATION PEDIATRIC  
10PCS

**OX20171A**



Rev.1 (15/01/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)