

# TEXAS 110

## Pulse oximeter



Portable device for non-invasive measuring of oxygen saturation of arterial blood, pulse rate and pulse strength.

### Specific features

The red light display allows easy reading  
 Judgment of exceeding the limits, Audio and visual alarm  
 Automatic shutdown when no pulse oxymetering for 3 minutes  
 Operates with 3 batteries AA. The operating time is up to 7 hours of continuous operation, with new batteries  
 Battery capacity indicator  
 Supplied with adult finger sensor

### Technical data

Dimensions	120x63x32 mm
SpO <sub>2</sub> measurement	35% - 99%
SpO <sub>2</sub> accuracy	70% - 99%, ± 2%; ≤69%, n.d.
SpO <sub>2</sub> Resolution	1%
Pulse measurement	30 - 250 BPM
Pulse error	30 - 250 BPM ± 2bpm
Pulse resolution	1 bpm
Weight	125g (without sensor and batt.)
Operative conditions	Temp. from 0 to 45°C Humid. from 15% to 95%
Storage conditions	Temp. from -20 a +60°C Humid. from 10% a 95%

Class IIa Med. Dev. according to Dir. 93/42/CEE  
 PULSEOXYMETER WITH SENSOR

PO00110

CND Classification Z1203020408

Distributed by Spencer Italia s.r.l.

Rev.1 (10/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)