

TEXAS 111

Compact pulse oximeter



Portable device for non-invasive measuring of oxygen saturation of arterial blood and pulse rate.

Specific features

- The display provides easy reading
- The operation is automatic and does not require the use of buttons
- Low battery visual warning
- Silicone case included

Technical data

| | |
|------------------------------|---|
| Dimensions | 57x34x31 mm |
| SpO ₂ measurement | 35% - 100% |
| SpO ₂ accuracy | 70% - 99%, ± 2%; ≤69%, n.d. |
| Pulse measurement | 30 - 250 BPM |
| Pulse error | 30 - 250 BPM ± 2% |
| Weight | 57g (batteries included) |
| Batteries | 2 type AAA |
| Applied part | Type B |
| Batteries operating time | About 30 hours |
| IP protection | IPX1 |
| Operative conditions | Temp. from 5 to 40°C Humid. from 15% to 90% |
| Storage conditions | Temp: from -20 to +70°C Humid. from 15% to 93% |

Compliant with

IEC 60601-1-2

Class IIa Med. Dev. according to Dir. 93/42/CEE

TEXAS 111 COMPACT
PULSEOXYMETER WITH
INTEGRATED SENSOR



PO00111C

CND Classification

Z1203020408

Distributed by Spencer Italia s.r.l.

Rev.5 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it