

PRX 04

Pressure reducer



Pressure reducer for direct connection to the oxygen tank. The pressure gauge allows easy reading of the pressure inside the oxygen tank.

Specific features

- Compact and reliable
- In put connection standard UNI
- Safety valve
- Output hose holder $\varnothing 7$ mm
- Integrated sintered filter

Technical data

Dimensions	100x135x40 \pm 5 mm
Weight	0,56 \pm 0,2 kg
Max input pressure	200 bar
Output pressure	Precalibrated max 4 bar
Output flow	6 l/min
Materials	Chromed brass, Steel, Bronze, Vyton
Manometer scale	0 – 300 bar precision class 2,5

Compliant with

ISO 10524-1:2006

Class IIb Med.Device according to Dir. 93/42/CEE

REDUCTION GEAR PRESSURE
PISTON B 3,5 6l/m.

PR00304B

CND Classification



Z120309

Registration number

84047

Distributed by Spencer Italia S.r.l

Rev.2 (15/01/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it