

## PRX 08

### Pressure reducer with flowmeter



Pressure reducer for direct connection to the oxygen tank. The pressure gauge allows easy reading of the pressure inside the oxygen tank.

#### Specific features

Input connection standard NF  
 Safety valve  
 Output with hose holder  $\varnothing 6\text{mm}$   
 Sintered bronze filter integrated in the device  
 Flowmeter with steps with ABS knob

#### Technical data

|                         |                                     |
|-------------------------|-------------------------------------|
| Dimensions              | 100x110x130 max $\pm 5$ mm          |
| Weight                  | 0,78 $\pm$ 0,2 kg                   |
| Max input pressure      | 200 bar                             |
| Output pressure         | Precalibrated max 4 bar             |
| Flow adjustment (l/min) | 0-1-1,5-2-3-4-6-9-12-15             |
| Materials               | Chromed brass, Steel, Bronze, Vytan |
| Manometer               | 0 – 300 bar precision class 2,5     |

#### Compliant with

ISO 10524-1:2006

Class IIb Med.Device according to Dir. 93/42/CEE

PRES.REG.DISP.15L/MIN NF  
 W/GAUGE GUARD OUTLET  
 $\varnothing 6\text{MM}$



**PR00316B**

**CND Classification** Z120309

**Registration number** 84055

Distributed by Spencer Italia S.r.l

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)