

Pressure reducer



Pressure reducer for direct connection to the oxygen tank. The pressure gauge allows easy reading of the pressure inside the oxygen tank.

Specific features

Input connection standard NF
Safety valve
Flowmeter with steps with ABS knob

Technical data

Dimensions	135x120x45 ± 10 mm
Weight	0,68 ± 0,2 kg
Max input pressure	200 bar
Output pressure	Precalibrated max 4 bar
Materials	Chromed brass, Steel, Bronze, Vytan
output flow	140 l/min
Output	¼" G
Manometer	0 – 300 bar precision class 2,5

Compliant with

ISO 10524-1:2006

Class IIb Med.Device according to Dir. 93/42/CEE

PRESSURE REDUCER 140 L/MIN
3,5 BAR NF STANDARD

PR00331

CND Classification



Z120309

Distributed by Spencer Italia S.r.l

Rev.10(24/05/2019)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it