

## Pressure reducer with flowmeter



Pressure reducer for direct connection to the oxygen tank. The pressure gauge allows easy reading of the pressure inside the oxygen tank.

### Specific features

- Input connection standard UNI
- Safety valve
- Output with hose holder  $\varnothing 6$ mm and UNI socket
- Sintered bronze filter integrated in the device
- Flowmeter with steps with ABS knob

### Technical data

Dimensions	140x110x130 max $\pm 10$ mm
Weight	1,02 $\pm 0,2$ kg
Max input pressure	200 bar
Output pressure	Precalibrated max 4 bar
Flow adjustment (l/min)	0-1-1,5-2-3-4-6-9-12-15
Materials	Chromed brass, Steel, Bronze, Vytan
Socket output flow	140 l/min
Manometer	0 – 300 bar precision class 2,5

### Compliant with

ISO 10524-1:2006

Class IIb Med.Device according to Dir. 93/42/CEE

PRESSURE REDUCER WITH  
DISPENSER AND OXYGEN  
SOCKET UNI



**PR00333B**

CND Classification

Z120309

Distributed by Spencer Italia S.r.l

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)