

PRX05

Pressure reducer

With trasduttore



Two-stage pressure reducer with fixed setting suitable for direct connection to the tank. The integrated transducer allows sending pressure readings to a monitoring device installed inside the medical vehicle.

Specific Characteristics

Compact and reliable
Inlet connection standard UNI
Safety valve
Outlet connection ¼" with ø7mm hose
Sintered bronze filter integrated in the device

Technical data

Dimensions	115x127x125 mm
Weight	1,42 kg
Max inlet pressure	230 bar
Output pressure	Adjustable from 0 to 7 bar
Flow	150 l/min
Filter	60 µ
Power supply	DC 8...30V
Output signal	0...5V, 3 fili
Trasduttore cable	48 cm
Materials	Chromed brass, Steel, Bronze, Vyton

Class IIb Med.Device according to Dir. 93/42/CEE

PRX05 WITH PRESSURE
TRASDUTTORE 300BAR
PR00336B



CND Classification	Z120309
Registration number	84047/R

Distributed by Spencer Italia S.r.l

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it