

Pressure reducer



Pressure reducer for direct connection to the oxygen tank. The pressure gauge allows easy reading of the pressure inside the oxygen tank.

Specific features

- Compact and reliable
- In put connection standard BS
- Safety valve
- Output hose holder $\varnothing 7\text{mm}$
- Integrated sintered filter

Technical data

Max input pressure	200 bar
Output pressure	Precalibrated max 4 bar
Output flow	140 l/min
Materials	Chromed brass, Steel, Bronze, Vytan
Manometer scale	0 – 300 bar precision class 2,5

Compliant with

ISO 10524-1:2006

Class IIb Med.Device according to Dir. 93/42/CEE

PRESS.RED.PIST.3,5-140L/m-
SHORT-W/COVER TUB.S.Ø7MM
PR00339



Distributed by Spencer Italia S.r.l

Rev.0 (28/02/2019)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it