

PRX 08

Pressure reducer with integrated flowmeter



Pressure reducer with fixed output pressure suitable for direct connection to the oxygen tank. The pressure gauge allow to read the remaining pressure inside the oxygen tank. The output flow can be adjusted thanks to the integrated flowmeter.

Specific Characteristics

Compact and reliable
 Input connection standard DIN
 Safety valve
 Output connection \varnothing 6mm hose
 Sintered bronze filter integrated in the device

Technical data

Dimensions	105x115x85 \pm 5 mm
Weight	720 \pm 20 g
Max inlet pressure	200 bar
Output pressure	Factory set to 4 bar
Available flow settings	0 – 1 – 1,5 – 2 – 3 – 4 – 6 – 9 – 12 – 15 l/min
Materials	Chromed brass, Steel, Bronze, Vytan
Manometer scale	0 – 300 bar

Compliant with

ISO 10524-1:2006

Class IIb Med.Device according to Dir. 93/42/CEE

PRESSURE REDUCER O2 DIN 15
 l/m-SHORT-HOSE \varnothing 6mm

PR00346B

CND Classification



Z120309

Distributed by Spencer Italia S.r.l

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it