

Pressure reducer with flowmeter



Pressure reducer for direct connection to the oxygen tank. The pressure gauge allows easy reading of the pressure inside the oxygen tank.

Specific features

Input connection standard NF
 Safety valve
 Output with hose holder $\varnothing 6$ mm and NF socket
 Sintered bronze filter integrated in the device
 Flowmeter with steps with ABS knob

Technical data

Dimensions	100x110x130 max ± 10 mm
Weight	0,85 $\pm 0,2$ kg
Max input pressure	200 bar
Output pressure	Precalibrated max 4 bar
Flow adjustment (l/min)	0-1-1,5-2-3-4-6-9-12-15
Materials	Chromed brass, Steel, Bronze, Vytan
Socket output flow	140 l/min
Manometer	0 – 300 bar precision class 2,5

Compliant with

ISO 10524-1:2006

Class IIb Med.Device according to Dir. 93/42/CEE

STEP PRESSURE REDUCER 0-15L/M NF WITH NF O2 OUTLET
PR00352



CND Classification

Z120309

Distributed by Spencer Italia S.r.l

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it