

PRX 07

Pressure reducer with fixed flow



Pressure reducer for direct connection to the oxygen tank. The pressure gauge allows easy reading of the pressure inside the oxygen tank.

Specific features

- Compact and reliable
- Input connection standard UNI
- Safety valve
- Output quick socket DIN
- Sintered bronze filter integrated in the device

Technical data

Dimensions	125x115x35 ± 5 mm
Weight	0,67 ± 0,2 kg
Max input pressure	200 bar
Output pressure	Precalibrated max 4 bar
Output flow	140 l/min ± 10%
Materials	Chromed brass, Steel, Bronze, Vytan
Pressure Gauge	0 – 300 bar precision class 2,5

Compliant with

ISO 10524-1:2006

Med.Device CE marked according to Dir. 93/42/CEE

RID.PRES.PIST.B3,5-140l/m-
CORT-C/COPR.PORT.Ø7mmUNI

PR00356



Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it