

## Pressure reducer with Click flowmeter and humidifier



Pressure reducer suitable for direct connection to the oxygen tank. The pressure gauge allow to read the remaining pressure inside the oxygen tank. The integration of the humidifier makes it more practical during oxygen therapy.

### Specific features

- Compact and reliable system
- Input connection standard DIN
- Safety valve
- The humidifier helps to reduce the dryness of mucosas during administration of oxygen
- Output hose holder  $\varnothing 6$ mm and DIN socket
- Sintered bronze filter integrated in the device

### Reference standard

- ISO 10524-1 – Pressure reducer
- UNI EN ISO 8185 - Humidifier
- DIN 477 nr 9 – Connection to oxygen tank

### Technical data

|                                   |   |
|-----------------------------------|---|
| Dimensions                        | 140x130x255 $\pm$ 20mm                          |
| Weight                            | 940 $\pm$ 30 g                                  |
| Max input pressure                | 200 bar   |
| Output pressure                   | Factory set max 4 bar                           |
| Flow adjustments                  | 0 – 1 – 1,5 – 2 – 3 – 4 – 6 – 9 – 12 – 15 l/min |
| Socket output flow                | 150 l/min                                       |
| Materials                         | Chromed brass, Steel, Bronze, Vyton, PC         |
| Manometer scale                   | 0 – 300 bar                                     |
| Humidifier max Volume             | 130 ml  |
| Humidifier max operating pressure | 500 kPa   |

Med. Dev. CE marked according to Dir. 93/42/CEE

PRESS.RED. WITH SELECTOR 0-  
15 L/M DIN W/HUM. AND O2  
DIN OUTLET

PR00360



Rev.1 (15/01/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)