

Pressure reducer with Click flowmeter



Pressure reducer suitable for direct connection to the oxygen tank. The pressure gauge allow to read the remaining pressure inside the oxygen tank.

Specific features

Compact and reliable system
 Input connection standard DIN
 Safety valve
 Output hose holder \varnothing 6mm or 9/16" thread
 DIN socket for supplementary outlet
 Sintered bronze filter integrated in the device

Technical data

Max input pressure	200 bar
Output pressure	Factory set max 4 bar
Flow adjustments	0 – 1 – 1,5 – 2 – 3 – 4 – 6 – 9 – 12 – 15 l/min
Socket output flow	150 l/min
Materials	Chromed brass, Steel, Bronze, Vyton, PC
Manometer scale	0 – 300 bar

Reference standard

ISO 10524-1 – Pressure reducer
 DIN 477 nr 9 – Connection to oxygen tank
 DIN 13260-2

Med. Dev. CE marked according to Dir. 93/42/CEE
 PRX08 DIN O2 W/DIN OULET
 W/HOSE NIPPLE
PR00365



Rev.0 (02/01/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it