

# KYOTO

## Cervical collar



Two – piece cervical collar, with hole for tracheostomized patient suitable for stabilization of cervical fractures.

### Specific Characteristics

Soft anti-decubitus material to increase patient comfort

Anatomical reinforcements of both valves to increase support

Strap adjustable closure

Holes to increase transpiration

Washable and reusable (for a single patient)

### Technical data

Size	S
Materials	Expanded PU, PE
Height <small>(Chin support - sternum)</small>	11 cm

Class I Med.Device according to Dir. 93/42/CEE

KYOTO - COLLAR SIZE S

QC11000B



CND Classification Y060312

Registration number 1215507

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)