

KYOTO

Cervical collar



Two – piece cervical collar, with hole for tracheostomized patient suitable for stabilization of cervical fractures.

Specific Characteristics

Soft anti-decubitus material to increase patient comfort

Anatomical reinforcements of both valves to increase support

Strap adjustable closure

Holes to increase transpiration

Washable and reusable (for a single patient)

Technical data

Size	M
Materials	Expanded PU, PE
Height <small>(Chin support - sternum)</small>	12 cm

Class I Med.Device according to Dir. 93/42/CEE

KYOTO - COLLAR SIZE M

QC11001B



CND Classification Y060312

Registration number 1215512

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it