

KYOTO

Cervical collar



Two – piece cervical collar, with hole for tracheostomized patient suitable for stabilization of cervical fractures.

Specific Characteristics

Soft anti-decubitus material to increase patient comfort

Anatomical reinforcements of both valves to increase support

Strap adjustable closure

Holes to increase transpiration

Washable and reusable (for a single patient)

Technical data

Size	XL
Materials	Expanded PU, PE
Height <small>(Chin support - sternum)</small>	13,5 cm

Class I Med.Device according to Dir. 93/42/CEE

KYOTO - COLLAR SIZE XL

QC11003B



CND Classification Y060312

Registration number 1238817

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it