

## PARIS

### One-piece soft cervical collar



One-piece cervical collar designed to minimize muscle contraction increasing neck stabilization and to relieve many neck diseases.

#### Specific features

- Light and comfortable
- Fast application thanks to the strap closure
- The cotton tissue increases skin transpiration
- The device contains latex
- Single patient use

#### Technical data

Size	L
Dimensions	590x130x25 ± 5mm
Material	PE exp., cotton
Weight	110 ± 5 g

#### Class I Med.Device according to Dir. 93/42/CEE

PARIS - SOFT CERVICAL COLLAR  
SIZE L

QC12002B

CND Classification Y060312

Registration number 1478640



Rev.0 (09/11/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)