

NEC YORK

One-piece cervical collar




One-piece cervical collar to be used in first aid and emergency situations. Suitable for the immobilization of the cervical spine in trauma patients.

Specific Characteristics

- Fast application thanks to the strap closure
- Access for tracheotomy
- Openings in the nuchal part for better transpiration
- Comfortable shaped padding

Technical data

Size	Infant
Materials	PE
Weight	85 g
Dimension A	125 mm
Dimension B	450 mm
Dimension C	80 mm



Accessories

- QC19010C - COLLAR BAG 1 – Bag for one-piece collars

Class I Med.Device according to Dir. 93/42/CEE

NEC YORK SIZE 1 (INFANT) -
1PIECE CERVIC.COLLAR



QC19001C

CND Classification

Y060312

Registration number

81010

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it