

NEC YORK

Kit 3 sizes of one-piece cervical collar



One-piece cervical collar to be used in first aid and emergency situations. Suitable for the immobilization of the cervical spine in trauma patients.

Specific Characteristics


- Fast application thanks to the strap closure
- Access for tracheotomy
- Openings in the nuchal part for better transpiration
- Comfortable shaped padding
- 3 sizes (pediatric, neckless, regular)
- Storage bag

Standard Accessories

- QC19010C - COLLAR BAG 1 – Bag for one-piece collars

Technical data

Sizes	2 (Pediatric) 3 (Neckless) 5 (Regular)
Materials	PE , PA
Kit weight (with bag)	430gr± 50gr
Dimension A (±5mm)	125 mm (Pediatric) 135 mm (Neckless) 145 mm(Regular)
Dimension B (±5mm)	445 mm(Pediatric) 555 mm (Neckless) 555 mm(Regular)
Dimension C (±5mm)	100 mm (Pediatric) 95 mm (Neckless) 115 mm(Regular)



Class I Med.Device according to Dir. 93/42/CEE

NEC YORK KIT 3SIZES(2,3,5)
WITH BAG 1PIECE C.COLL.

QC19013C

CND Classification



Y060312

Registration number

195214

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it