

NEC YORK

Kit 6 one-piece cervical collars



Nec York cervical collars are devices to be used to help maintain neutral alignment of the cervical spine.

Specific features


- Fast application thanks to the strap closure
- Access for tracheotomy
- Openings in the nuchal part for better transpiration
- Comfortable shaped padding
- 5 sizes (Infant, pediatric, neckless, 2x short, regular)
- Storage bag

Standard equipment

QC19010C - COLLAR BAG 1 – Bag for one-piece collars

Technical data

Materials	PE, PA
Kit weight (with bag)	760 ± 50 g
Dimension A (± 5mm)	125 mm (Infant)
	125 mm (Paediatric)
	135 mm (Neckless)
	130 mm (Short)
	145 mm (Regular)
Dimension B (± 5mm)	450 mm (Infant)
	445 mm (Paediatric)
	555 mm (Neckless)
	555 mm (Short)
	555 mm (Regular)
Dimension C (±5 mm)	80 mm (Infant)
	100 mm (Paediatric)
	95 mm (Neckless)
	80 mm (Short)
	115 mm (Regular)



Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it

QC19015C

NEC YORK KIT 6SIZES WITH BAG-1-PIECE CERV.COLL.

SET consisting of the following MD compliant with UE Reg. 2017/745

QC19001C - NEC YORK SIZE 1 (INFANT) - 1PIECE CERVIC.COLLAR	Reg. n. 81010	CND Y060312
QC19002C – NEC YORK SZ.2 - 1 PIECE RIGID CERVICAL COLLAR PEDI	Reg. n. 83084	CND Y060312
QC19020C – NEC YORK SIZE 3 (NECKLESS) - 1PIECE CERV. COLLAR	Reg. n. 83086	CND Y060312
QC19021C – NEC YORK SIZE 4 (SHORT) - ONE- PIECE CERV. COLLAR	Reg. n. 83088	CND Y060312
QC19022C – NEC YORK SIZE 5 (REGULAR) - 1PIECE CERV. COLLAR	Reg. n. 83090	CND Y060312

Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it