

NEC YORK

Kit 6 pezzi collari cervicali monovalva



I collari cervicali Nec York sono dispositivi da utilizzarsi per contribuire al mantenimento dell'allineamento neutro del rachide cervicale.

Caratteristiche specifiche


- Rapida applicazione grazie alla chiusura strap
- Accesso per tracheotomia
- Aperture nella parte nucale per migliorare l'aerazione
- Imbottitura sagomata confortevole
- 5 Misure (infant, paediatric, neckless, 2x short, regular)
- Sacca di trasporto

Accessori opzionali

- QC19010C - COLLAR BAG 1 – Sacca per collari monovalva

Dati tecnici

Materiali	PE, PA
Peso kit (con sacca)	760 ± 50 g
Dimensione A (± 5 mm)	125 mm (Infant)
	125 mm (Paediatric)
	135 mm (Neckless)
	130 mm (Short)
	145 mm (Regular)
Dimensione B (± 5 mm)	450 mm (Infant)
	445 mm (Paediatric)
	555 mm (Neckless)
	555 mm (Short)
	555 mm (Regular)
Dimensione C (± 5 mm)	80 mm (Infant)
	100 mm (Paediatric)
	95 mm (Neckless)
	80 mm (Short)
	115 mm (Regular)



Rev.0 (31/03/2021)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it

QC19015C

 NEC YORK - CERVICAL COLLARS KIT 6 PZ. (5 MISURE)
 C/SACCA

SET costituito dai seguenti DM conformi al Reg.2017/745

QC19001C - NEC YORK - COLLARE CERVICALE MIS.1 (INFANT)	Id. Reg. 81010	CND Y060312
QC19002C - NEC YORK - COLLARE CERVICALE MIS.2 (PAEDIATRIC)	Id. Reg. 83084	CND Y060312
QC19020C - NEC YORK - COLLARE CERVICALE MIS.3 (NECKLESS)	Id. Reg. 83086	CND Y060312
QC19021C - NEC YORK - COLLARE CERVICALE MIS.4 (SHORT)	Id. Reg. 83088	CND Y060312
QC19022C - NEC YORK - COLLARE CERVICALE MIS.5 (REGULAR)	Id. Reg. 83090	CND Y060312

Rev.0 (31/03/2021)

 COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it