

NEC YORK

One-piece cervical collar




Nec York cervical collars are devices to be used to help maintain neutral alignment of the cervical spine.

Specific features

- Fast application thanks to the strap closure
- Access for tracheotomy
- Openings in the nuchal part for better transpiration
- Comfortable shaped padding

Technical data

Size	Regular
Materials	PE
Weight	130 g
Dimension A	145 mm
Dimension B	555 mm
Dimension C	120 mm



Optional accessories

- QC19010C - COLLAR BAG 1 – Bag for one-piece collars

Class I MD compliant with UE Reg. 2017/745

NEC YORK SIZE 5 (REGULAR)
- 1PIECE CERV. COLLAR

QC19022C

CND Classification Y060312

Registration number 83090



Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it