

## NECK YORK

One-piece cervical collar




One-piece cervical collar to be used in first aid and emergency situations. Suitable for the immobilization of the cervical spine in trauma patients.

### Specific Characteristics

- Fast application thanks to the strap closure
- Access for tracheotomy
- Openings in the nuchal part for better transpiration
- Comfortable shaped padding

### Technical data

Size	Large
Materials	PE
Weight	120 g
Dimension A	135 mm
Dimension B	555 mm
Dimension C	130 mm



### Accessories

QC19010C - COLLAR BAG 1 – Bag for one-piece collars

Class I Med.Device according to Dir. 93/42/CEE

NEC YORK SIZE 6 - 1PIECE  
CERV. COLLAR (LARGE)

QC19024C

CND Classification



Y060312

Rev.1 (07/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)