

## NEC YORK

### Kit 6 one-piece cervical collars




Nec York cervical collars are devices to be used to help maintain neutral alignment of the cervical spine.

#### Specific features

- Fast application thanks to the strap closure
- Access for tracheotomy
- Openings in the nuchal part for better transpiration
- Comfortable shaped padding
- The kit includes 4 regular size and 2 paediatric size collars
- Storage bag

#### Technical data

Materials	PE, PA
Dimension A (±5mm)	125 mm (Paediatric) 145 mm (Regular)
Dimension B (±5mm)	445 mm (Paediatric) 555 mm (Regular)
Dimension C (±5mm)	100 mm (Paediatric) 115 mm (Regular)



#### Standard equipment

QC19010C - COLLAR BAG 1 – Bag for one-piece collars

#### QC19025C

NEC YORK KIT 6 COLLARS WITH BAG

SET consisting of the following MD compliant with UE Reg. 2017/745

QC19002C – NEC YORK SZ.2 - 1 PIECE RIGID CERVICAL COLLAR PEDI	Reg. n. 83084	CND Y060312
QC19022C – NEC YORK SIZE 5 (REGULAR) - 1PIECE CERV. COLLAR	Reg. n. 83090	CND Y060312

Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)