

NEC YORK

Kit 6 one-piece cervical collar



One-piece cervical collar to be used in first aid and emergency situations. Suitable for the immobilization of the cervical spine in trauma patients.

Specific features


- Fast application thanks to the strap closure
- Access for tracheotomy
- Openings in the nuchal part for better transpiration
- Comfortable shaped padding
- The kit includes 4 Regular size and 2 pediatric size collars
- Storage bag

Technical data

Materials	PE , PA
Dimension A (±5mm)	125 mm (Pediatric)
	145 mm (Regular)
Dimension B (±5mm)	445 mm (Pediatric)
	555 mm (Regular)
Dimension C (±5mm)	100 mm (Pediatric)
	115 mm (Regular)

Standard Accessories

QC19010C - COLLAR BAG 1 – Bag for one-piece collars



Class I Med.Device according to Dir. 93/42/CEE

NEC YORK KIT 6 COLLARS WITH BAG

QC19025C

CND Classification



Y060312

Rev.0 (02/09/2019)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it