


NEC UP

Multi-size adult collar



Nec Up cervical collars are devices to be used to help maintain neutral alignment of the cervical spine.

Specific features

- Quick application thanks to the strap closure
- Access for tracheotomy
- Padding to increase comfort
- Rear hole for ventilation and to allow inspections
- Chin support has 4 heights in order to make the device suitable for a wide range of cervical conformations
- The selected size can be easily identified thanks to the red colored marking visible on the front side
- After the size has been selected, the chin support can be easily fixed by pressing the two locking tabs
- Compatible with RX, MRI, TAC ^(Note)

Note: The compatibility is expressed in reference to the use of non-magnetic materials

Technical data

Obtainable chin support heights	Tall – 135 mm
	Regular – 125 mm
	Short – 115 mm
	Low – 105 mm
Circumference	575 mm
Material	PE
Weight	120 g

Class I MD compliant with UE Reg. 2017/745

NEC UP ADULT - MULTI-SIZE
COLLAR

QC19100B

CND Classification



Y060312

Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it