

## NEC UP

### Multi-size pediatric collar



One piece multi-size cervical collar designed for use in first aid and emergency situations. It's identified as a solution for the immobilization of the cervical spine in trauma patients.

#### Specific Characteristics

- Quick application thanks to the strap closure
- Access for tracheotomy
- Padding to increase comfort
- Rear hole for ventilation and to allow inspections
- Chin support has 3 heights in order to make the device suitable for a wide range of cervical conformations
- The selected size can be easily identified thanks to the red colored marking visible on the front side
- After the size has been selected, the chin support can be easily fixed by pressing the two locking tabs
- Compatible with RX, MRI, TAC <sup>(Note)</sup>

Note: The compatibility is expressed in reference to the use of non-magnetic materials

#### Technical data

Obtainable chin support heights	Petite – 80 mm
	Pedi – 70 mm
	Infant – 60 mm
Circumference	465 mm
Material	PE
Weight	85 g

Class I Med.Device according to Dir. 93/42/CEE

NEC UP PAEDIATRIC - MULTI-SIZE  
COLLAR

QC19101B

CND Classification



Y060312

Rev.4 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)