

DALLAS

One-piece cervical collar




One-piece cervical collar to be used in first aid and emergency situations for a single patient. They are suitable for the immobilization of the cervical spine in trauma patients.

Specific Characteristics

- Light and comfortable
- Fast application thanks to the strap closure
- The white color allows an immediate cleanliness visual feedback

Technical data

Dimension A	120 mm
Dimension B	95 mm
Dimension C	405 mm
Dimension D	65 mm
Material	PE exp.
Weight	16 g



Class I Med. Device according to Dir. 93/42/CEE

DALLAS SIZE XS -SEMIRIGID ONE-PIECE CERV. COLLAR

QC22301A

CND Classification



Y060312

Registration number

92291

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it