

DALLAS

Semi-rigid one-piece cervical collar




Dallas cervical collars are devices to be used to help maintain neutral alignment of the cervical spine.

Specific features

- Light and comfortable
- Fast application thanks to the strap closure
- The white color allows an immediate cleanliness visual feedback

Technical data

Dimension A	120 mm
Dimension B	95 mm
Dimension C	405 mm
Dimension D	65 mm
Material	Expanded PE
Weight	16 g



Class I MD compliant with UE Reg. 2017/745

DALLAS SIZE XS -SEMIRIGID ONE-PIECE CERV. COLLAR

QC22301A

CND Classification



Y060312

Registration number

92291

Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it