

DALLAS

Kit 5 semi-rigid one-piece cervical collars



Dallas cervical collars are devices to be used to help maintain neutral alignment of the cervical spine.

Specific features

- Light and comfortable
- Fast application thanks to the strap closure
- The white color allows an immediate cleanliness visual feedback

Optional accessories

- QC19010C - COLLAR BAG 1 – Bag for one-piece collars

Technical data

Included sizes	XS, S, M, L, XL
Materials	Expanded PE
Weight of the kit	200 g

QC22370A

DALLAS KIT FIVE SIZES-SEMI-RIGID 1PIECE CERV.COLL.

SET consisting of the following MD compliant with UE Reg. 2017/745

QC22301A – DALLAS SIZE XS - SEMI-RIGID ONE-PIECE CERV. COLLAR	Reg. n. 92291	CND Y060312
QC22302A – DALLAS SIZE S - SEMI-RIGID ONE-PIECE CERV. COLLAR	Reg. n. 92293	CND Y060312
QC22303A - DALLAS SIZE M - SEMI-RIGID ONE-PIECE CERV. COLLAR	Reg. n. 92298	CND Y060312
QC22304A – DALLAS SIZE L - SEMI-RIGID ONE-PIECE CERV. COLLAR	Reg. n. 92300	CND Y060312
QC22305A – DALLAS SIZE XL - SEMI-RIGID ONE-PIECE CERV. COLLAR	Reg. n. 92303	CND Y060312

Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it