


CERVI FLEX

One-piece cervical collar



One-piece cervical collar to be used in first aid and emergency situations. They are suitable for the immobilization of the cervical spine in trauma patients.


Specific Characteristics

Light and comfortable
Fast application thanks to the strap closure
Extremely easy to clean

Technical data

Size A	105 mm
Size B	590 mm
Size C	70 mm
Materials	Vinyl coated PU
Weight	85 ± 5 g

Dimensions subjected to tolerances of ± 5mm



Class I Med.Device according to Dir. 93/42/CEE

CERVI FLEX SIZE M - ONE-PIECE
CERVICAL COLLAR

QC22503A

CND Classification Y060312

Registration number 100545



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it