

CERVI FLEX

One-piece cervical collar



One-piece cervical collar to be used in first aid and emergency situations. They are suitable for the immobilization of the cervical spine in trauma patients.


Specific Characteristics

- Light and comfortable
- Fast application thanks to the strap closure
- Extremely easy to clean

Technical data

Size A	80 mm
Size B	460 mm
Size C	55 mm
Materials	Vinyl coated PU
Weight	60 ± 5 g

Dimensions subjected to tolerances of ± 5mm



Class I Med.Device according to Dir. 93/42/CEE

CERVI FLEX SIZE XS - ONE-PIECE
CERVICAL COLLAR

QC22504A

CND Classification



Y060312

Registration number

99955

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it