

CERVI FLEX

One-piece cervical collar



One-piece cervical collar to be used in first aid and emergency situations. They are suitable for the immobilization of the cervical spine in trauma patients.

Specific Characteristics

- Light and comfortable
- Fast application thanks to the strap closure
- Extremely easy to clean

Technical data

Size A	125 mm
Size B	655 mm
Size C	95 mm
Materials	Vinyl coated PU
Weight	105 ± 5 g

Dimensions subjected to tolerances of ± 5mm



Class I Med.Device according to Dir. 93/42/CEE

CERVI FLEX SIZE XL - ONE-PIECE
CERVICAL COLLAR
QC22505A



CND Classification Y060312
Registration number 100548

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it