


CERVI FLEX

Kit 5 one-piece semirigid cervical collars



One-piece cervical collar to be used in first aid and emergency situations. They are suitable for the immobilization of the cervical spine in trauma patients.

Specific Characteristics

Light and comfortable
Fast application thanks to the strap closure
Extremely easy to clean

Technical data

Included sizes	XS, S, M, L, XL
Materials	Vinyl coated PU
Weight of the kit	420 ± 10 g

Equipment

QC19010C - COLLAR BAG 1 – Bag for one piece collars

Class I Med.Device according to Dir. 93/42/CEE

CERVI FLEX KIT FIVE SIZES - ONE-PIECE CERV. COLLAR

QC22506A

CND Classification Y060312

Registration number 100551



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it