

YES

Cervical collar



YES cervical collars are devices to be used to help maintain neutral alignment of the cervical spine.

Specific features

- Anatomic shape of the occipital and mandibular zones
- Fast application thanks to the strap closure
- Access for tracheotomy
- Dedicated opening for the control of carotid pulse
- Comfortable shaped padding
- The foamed construction material confers lightness

Technical data

Size	S
Materials	Expanded PE, nylon

Class I MD compliant with UE Reg. 2017/745

YES - TWO PIECES CERVICAL
COLLAR SIZE S

QC30101A

CND Classification



Y060312

Registration number

56293981

Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it