

YES

Cervical collar



YES cervical collars are devices to be used to help maintain neutral alignment of the cervical spine.

Specific features

Anatomic shape of the occipital and mandibular zones

Fast application thanks to the strap closure

Access for tracheotomy

Dedicated opening for the control of carotid pulse

Comfortable shaped padding

The foamed construction material confers lightness

Technical data

Size	L
Materials	Expanded PE, nylon

Class I MD compliant with UE Reg. 2017/745

YES - TWO PIECES CERVICAL
COLLAR SIZE L

QC30103A

CND Classification Y060312

Registration number 565931



Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it