

YES

Cervical collars - Kit 4 sizes with bag



YES cervical collars are devices to be used to help maintain neutral alignment of the cervical spine.

Specific features

Anatomic shape of the occipital and mandibular zones

Fast application thanks to the strap closure

Access for tracheotomy

Dedicated opening for the control of carotid pulse

Comfortable shaped padding

The foamed construction material confers lightness

The kit includes a carry bag with zip closure and possibility to be hung

Technical data

Included sizes	XS, S, M, L
Materials	Expanded PE, nylon
Weight of the kit	545 g
Bag dimensions	420 x 330 mm

QC30104A

YES - CERVICAL COLLARS KIT 4 SIZES WITH BAG

SET consisting of the following MD compliant with UE Reg. 2017/745

QC30100A – YES - TWO PIECES CERVICAL COLLAR SIZE XS	Reg. n. 565927	CND Y060312
QC30101A – YES - TWO PIECES CERVICAL COLLAR SIZE S	Reg. n. 565929	CND Y060312
QC30102A – YES - TWO PIECES CERVICAL COLLAR SIZE M	Reg. n. 565930	CND Y060312
QC30103A – YES - TWO PIECES CERVICAL COLLAR SIZE L	Reg. n. 565931	CND Y060312

Rev.0 (31/03/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it