

## FIXO SPLINT

### Rigid splint



Fixo Splint splints are devices used to immobilise limbs in the event of suspected injuries or fractures of the limbs, to avoid additional strain during the pre-hospital phase.

#### Specific features

- Low storage space required
- Quick application thanks to the strap belt system
- Heat sealed coating to ensure easy sanitation
- Soft padding to increase patient comfort

#### Technical data

Dimensions <sup>(1)</sup>	120 x 220 x 310 ± 10 mm
Materials	PVC, wood
Weight	410 g

<sup>(1)</sup> Width x Depth x Height (as shown in the picture)

Class I MD compliant with UE Reg. 2017/745

FIXO SPLINT SIZE XS - RIGID  
SPLINTS

**QC70010A**

**CND Classification**



V0804

**Registration number**

103665

Rev.0 (04/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)