

FIXO SPLINT KIT

Rigid splints – Kit 2 sizes with bag



Fixo Splint splints are devices used to immobilise limbs in the event of suspected injuries or fractures of the limbs, to avoid additional strain during the pre-hospital phase.

Specific features

- Low storage space required
- Quick application thanks to the strap belt system
- Heat sealed coating to ensure easy sanitation
- Soft padding to increase patient comfort

Standard equipment

- Orange PVC carry bag

Technical data

Splint	Dimensions ⁽¹⁾ (mm)	Weight (g)
Size S	120 x 230 x 500 ± 10	640
Size L	150 x 230 x 730 ± 10	870
Materials	PVC, wood	
Weight of the kit	2 kg	
Bag dimensions	800 x 430 mm	
Bag weight	450 g	

⁽¹⁾ Width x Depth x Height (as shown in the picture)

QC70090A

FIXO SPLINT KIT TWO SIZES - RIGID SPLINTS

SET consisting of the following MD compliant with UE Reg. 2017/745

QC70020A – FIXO SPLINT SIZE S - RIGID SPLINTS	Reg. n. 103660	CND V0804
QC70050A – FIXO SPLINT SIZE L - RIGID SPLINTS	Reg. n. 103657	CND V0804

Rev.0 (04/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it