

## CDK Kit

### Tapered vacuum mattress with pump and bag



CDK vacuum mattresses are devices primarily intended to ensure immobilization of the whole body during patient transport.

#### Specific features

Robust, non-absorbent and self-extinguishing fabric  
 Perimetral rope with soft handles for easy transport  
 Tapered shape to improve adherence to the patient  
 9 carrying handles make it easier to transport and allow better adaptation in various rescue scenarios  
 Comes with ABS double effect pump and carry bag

#### Accessories

ST00593A	K Belt 3 – 3 belts 1 piece with bag
QM22199A	Repair kit

#### Technical data

Length	2150 ± 30 mm
Width	900 ± 30 mm
Number of handles	9
Materials	PVC, Nylon, Brass, Polystyrene
Weight (not including pump)	4,1 ± 0,3 kg
Load capacity	150 kg

Class I MD compliant with UE Reg. 2017/745

MOD.S018 KIT QM22300 +  
 QM22109 + QM22170



QM22092A

CND Classification V0804

Registration number 794749

Rev.0 (27/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)