

Res-Q-Matt Plus

Vacuum mattress



Res-Q-Matt Plus vacuum mattresses are devices primarily intended to ensure immobilization of the whole body during patient transport.

Specific features

Extremely resistant fabric
Integrated handles allow a firm grip

Technical data

Length	2140 ± 20 mm
Width	940 ± 10 mm
Filling material	Polystyrene
Number of handles	8
Materials	PVC, Nylon, Brass, Polystyrene
Weight (not including pump)	5 ± 0,5 kg
Load capacity	150 kg

Accessories

QM22120A	Aluminum pump
QM22125A	Aluminum double effect pump
QM22170A	Carry bag for vacuum mattresses
QM22199A	Repair kit

Class I MD compliant with UE Reg. 2017/745

RES-Q-MATT PLUS POLYSTYRENE-
PROFESS. VAC. MATTR.



QM22101A

CND Classification

V0804

Rev.0 (27/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it