

Res-Q-Matt

Compact vacuum mattress



Vacuum mattresses are devices for the immobilization of trauma patients. The vacuum system allows different levels of immobilization. The more rigid helps the right positioning of the patient; softer ones allow to mitigate the low-frequency vibrations typically found on emergency vehicles during transport.

Specific features

Extremely resistant fabric
Removable cover to make the cleaning procedures easier

Technical data

Length	2060 ± 20mm
Width	800 ± 10mm
N° of handles	8
Materials	PVC, Nylon
Weight	7,8 ± 0,5 kg

Included equipment

QM22120A	Aluminum Pump
QM22170A	Carry bag for vacuum mattresses

Class I Med.Device according to Dir. 93/42/CEE

RES-Q-MATT - ORANGE VAC. MAT.
W/BAG AND QMX120

QM22206A

CND Classification



V0804

Rev.0 (25/02/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it