

QMX 130

Aluminum vacuum pump



Compact vacuum pump, ideal for use with vacuum splints.

Specific features

The aluminum chamber makes the device strong and light
 Characterized by small dimensions, it can be easily stored
 in any vacuum splints kit
 Rounded knob for firm grip

Technical data

Dimensions	Ø55x320 mm
Tube internal diameter	8 mm
Tube length	255 ± 10mm
Materials	Al, PVC
Weight	400 ± 10 g

QMX 130 - ALUMINIUM PUMP FOR
 VACUUM SPLINTS
QM22130A

Rev.2 (31/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it