

Res-Q-Matt

Compact vacuum mattress



Res-Q-Matt vacuum mattresses are devices primarily intended to ensure immobilization of the whole body during patient transport.

Specific features

Extremely resistant fabric
Removable cover to make the cleaning procedures easier

Technical data

Length	2060 ± 20 mm
Width	800 ± 10 mm
Number of handles	8
Materials	PVC, Nylon, Polystyrene
Weight (not including pump)	7,8 ± 0,5 kg
Load capacity	150 kg

Accessories

QM22109A	ABS double effect pump
QM22120A	Aluminum pump
QM22170A	Carry bag for vacuum mattresses
QM22199A	Repair kit
ST00593A	K BELT 3 – Three belts STX 597

Class I MD compliant with UE Reg. 2017/745

RES-Q-MATT ORANGE - COMPACT VACUUM MATTRESS	
QM22206A	
CND Classification	V0804
Registration number	102652

Rev.0 (27/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it